MED EL

Fixation Clip

Ms010103

REF 09917

hearLIFE

AW30084_3.0

English 2
Deutsch 4
Français 6
Español 8
Português 10
Italiano 12
Nederlands 14
Svenska 16
Polski 18
Norsk 20
Suomi 22
Ελληνικα 24
Türkçe 26
Dansk 28
Magyar 30
Română 32
Slovenčina 34
Čeština 36
Български Език 38
Eesti Keel 40
Latviešu Valoda 42
Lietuvių Kalba 44
Slovenščina 46
Hrvatski 48
Русский 50
Srpski Jezik 52
한국어 54
中文 56
中文 58
61 עברית
63 العربية

Instructions for use


Fig. 1 Fixation Clip


Fig. 2 Distal part of the electrode

Fig. 3 Placement of the Fixation Clip on the bony bridge

Fig. 4 Bone fixation clip crimped on the bony bridge

Fig. 5 Electrode fixation clip closed over the electrode

Description of the device

The Fixation Clip is a single-use long-term passive implantable device to be implanted during MED-EL Cochlear Implant surgery. The Fixation Clip has two openings. The larger opening is the bone fixation clip and the smaller opening is the electrode fixation clip (see Fig. 1).

Intended use

The Fixation Clip is intended to be implanted by ENT surgeons during cochlear implant surgery to stabilise the electrode lead inside the middle ear cavity. The Fixation Clip should be used during cochlear implant surgery after the implant is immobilised and the electrode was inserted into the cochlea.
The Fixation Clip is intended for single use only.

Indication

The Fixation Clip is an implantable device and may be implanted in patients during MED-EL Cochlear Implant surgeries with all MED-EL electrodes with a diameter of 1.3 mm at the distal part of the electrode lead (see Fig. 2).

Contraindications

The Fixation Clip shall not be used with patients known to be intolerant of titanium.
The Fixation Clip shall not be used in the case that the bony bridge is not suitable for placing the Fixation Clip.

Materials

Titanium ASTM F67 Gr.2

Instructions for use

1. Before inserting the electrode array a bony bridge between the anterior tympanotomy and the posterior tympanotomy (facial recess) shall be created. The thickness of the bony bridge shall be created in such a way that it fits to the size of the bone fixation clip (1-2 mm is recommended) (see Fig. 1).
2. The Fixation Clip shall be implanted after the implant is immobilised and the electrode was inserted into the cochlea.
3. The openings are not of the same size. The larger opening of the Fixation Clip is the bone fixation clip and the smaller opening of the Fixation Clip is the electrode fixation clip (see Fig. 1).
4. First the larger clip (bone fixation clip) shall be placed on the bony bridge (see Fig. 3).
5. Next, the bone fixation clip shall be properly crimped to the bony bridge (see Fig. 4).
6. After fixing the bone fixation clip, the electrode fixation clip shall be closed over the electrode lead so that both ends of the electrode fixation clip come together and the clip is completely closed (see Fig. 5).
 - Care should be taken that during the fixation of the bone fixation clip the Incudostapedial joint is not harmed.
 - Care should be taken that following fixation of the bone fixation clip, the movements of the ossicles are not inhibited.

Undesirable side effects
Risks related to surgery

Possible post-operative side effects include the following: Dizziness, increased vertigo, delay in scar healing, an impaired sense of taste, potential for swallowing difficulties, numbness, increased tinnitus, stimulation of the facial nerve, temporary pain and uncomfortable sounds during stimulation.

Sterility

The Fixation Clip has been subjected to a validated gamma irradiation sterilisation process and is supplied in sterile packaging. Once the sterile packaging has been opened, the Fixation Clip cannot be resterilised. Do not use if the sterile packaging is damaged. The Fixation Clip is for single use only. Do not remove from sterile packaging until required.

Storage, shipment and disposal

The Fixation Clip shall only be shipped between –20 °C and +55 °C and stored in its sterile packaging at room temperature protected from moisture and direct heat.
Each device must be implanted before the use-by date specified on the package.
Dispose in accordance with federal, state and local regulations.

Information about use
General precautions and warnings

- The Fixation Clip shall not be altered and shall only be used as intended.
- The Fixation Clip shall be used exclusively with MED-EL electrodes that have a diameter of 1.3 mm at the distal part of the electrode lead.
- The distal part of the electrode lead is the part between the marker ring and the stimulator housing (see Fig. 2).
- Peel-off labels are included on the product packaging of the Fixation Clip. These labels are intended to be placed on the MED-EL Cochlear Implant Registration Card and in the patient's record.
- The Fixation Clip is a very small medical device and should be handled with appropriate care.
- Sterility of the Fixation Clip must be ensured at all times.
- The Fixation Clip shall not be used in cases of opened or damaged packaging.
- The Fixation Clip shall not be used after the expiry date.

Surgical precautions and warnings
Risks related to surgery

- Cochlear implant surgery is comparable to middle ear surgery with additional access to the inner ear. The normal risks of surgery and general anaesthesia are applicable. Primary surgical risks include the following: infection, inflammation, necrosis, haematoma, leakage of CSF, damage to the facial nerve, pain, scarring of the wound, swallowing difficulties and complications related to general anaesthesia. Additionally, meningitis can be a rare post-operative complication, but has the potential to be serious. The risk of meningitis may be reduced, for example by vaccination, antibiotic cover, and surgical technique.
- If available, facial nerve monitoring is recommended and if carried out, neural muscular blockade should be avoided.
- Accidental bending of the Fixation Clip during removal from its packaging must be avoided in order to prevent functional damage.
- Avoid excessive pressure on auditory structures of the middle ear.
- When implanting the Fixation Clip care must be taken that first the bone fixation clip is affixed.
- Care should be taken not to squeeze or damage the electrode.

Warranty

The warranty depends on the legal requirement. The warranty period begins from the date of delivery.

MRI caution

The Fixation Clip is intended to be used with a MED-EL Cochlear Implant. Therefore, the precautions and warnings for the implant apply to the Fixation Clip. Please refer to the MRI information provided in the implant instructions for use.

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MED-EL Elektromedizinische Geräte GmbH